



PRODUCT SPECIFICATION

α -Fetoprotein (AFP), Cell Culture

Catalog Number: A2414 **Part Number:** 90085

Source: Human Hepatoma Cell Line

Form

White powder, lyophilized from 10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4 \pm 0.1

Purity Specification

\geq 95% by SDS-PAGE

Protein / Content Specification

Report protein concentration by immunoassay prior to lyophilization (indicate method)

Characteristics Specification

Indicate calibration standard for immunoassay, if available

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.