



## PRODUCT SPECIFICATION

### C-Reactive Protein (CRP), Plasma, $\geq 80\%$

**Catalog Number:** C0123      **Part Number:** 90093

**Source:** Human Plasma

**Form**

Clear to hazy liquid - Buffer: 20mM Tris, 280 mM Sodium Chloride, 5mM Calcium Chloride, 0.1% Sodium Azide, pH 8.0  $\pm$  0.1

**Purity Specification**

$\geq 80\%$  by SDS-PAGE or equivalent

**Protein / Content Specification**

$\geq 2.0$  mg/ml by nephelometry

**Storage / Shipping**

Product is shipped using a refrigerated gel pack (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*