



PRODUCT SPECIFICATION

C-Reactive Protein (CRP), Pleural Fluid, $\geq 95\%$

Catalog Number: C0128 **Part Number:** 90429

Source: Human Pleural Fluid

Form

Clear to hazy liquid - Buffer: 20mM Tris, 280 mM Sodium Chloride, 5mM Calcium Chloride, 0.1% Sodium Azide, pH 8.0 \pm 0.1

Purity Specification

$\geq 95\%$ by SDS-PAGE

Protein / Content Specification

≥ 2.0 mg/ml (specify method)

Storage / Shipping

Product is shipped using a refrigerated gel pack (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.