



PRODUCT SPECIFICATION

Folate Binding Protein

Catalog Number: F0524 **Part Number:** 90044

Source: Bovine Milk

Form

Clear liquid - Buffer: 120mM Sodium Phosphate, 30mM Sodium Chloride, pH 7.4 ± 0.1

Purity Specification

Report results and method

Protein / Content Specification

Report protein concentration. This is determined indirectly using the measured activity in µg/ml and the binding capacity of 8 µg folic acid/mg folate binding protein

Activity Specification

Report folate binding ability in µg/ml

One mg of active folate binding protein will bind 8µg of folic acid

Storage / Shipping

Product is shipped using dry ice (***not temperature monitored***). On arrival, store frozen below -20°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.