



PRODUCT SPECIFICATION

Follicle Stimulating Hormone (FSH), Human, α -Subunit

Catalog Number: F0714 **Part Number:** 90106

Source: Human Pituitary Glands

Form

Lyophilized powder

Purity Specification

≥ 98% by SDS-PAGE

Characteristics Specification

This product is bulk processed and packaged by dry weight

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.