



PRODUCT SPECIFICATION

Myoglobin, Human Heart, $\geq 95\%$

Catalog Number: M0724 **Part Number:** 90026

Source: Human Heart

Form

Clear red-brown liquid - Buffer: 20mM Tris, 1mM EDTA, 50% Glycerol, pH 8.5 \pm 0.1

Purity Specification

$\geq 95\%$ by SDS-PAGE

Protein / Content Specification

≥ 2.0 mg/ml by Lowry Protein Assay (BSA standard)

Storage / Shipping

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen between -10°C to -25°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.