



PRODUCT SPECIFICATION

Prostatic Acid Phosphatase, Human, Liquid

Catalog Number: P0524 **Part Number:** 90219

Source: Human Seminal Fluid

Form

Clear liquid - Buffer: 10mM Tris, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 ± 0.1

Purity Specification

≥ 98% by SDS-PAGE

Protein / Content Specification

≥ 1.0 mg/ml by Lowry Protein Assay, BSA standard

Storage / Shipping

Product is shipped using a refrigerated gel pack (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.