



## **PRODUCT SPECIFICATION**

### **Plate Coating Stabilizer, no preservative**

**Catalog Number:** P3323      **Part Number:** 90227

#### **Form**

Clear golden liquid containing a proprietary mixture of non-toxic components, 0.2µm filtered. This product does not contain preservatives and must be handled under sterile conditions.

#### **Characteristics Specification**

Plate Coating Stabilizer is a proprietary mixture of components designed to increase the stability of immobilized antigens and antibodies. This solution is supplied ready-to-use and in many cases can replace current blocking solutions. The sterile liquid is supplied without preservatives and must be handled using aseptic techniques. Alternatively, a preservative compatible with the assay may be added.

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*