



PRODUCT SPECIFICATION

Sex Hormone Binding Globulin (SHBG), High Purity

Catalog Number: S1725 **Part Number:** 90242

Source: Human Serum

Form

Clear liquid - Buffer: 20mM Tris-HCl, 10% Glycerol, pH 7.4 ± 0.1

Purity Specification

≥ 95% by SDS-PAGE

Report contaminants

Protein / Content Specification

Report protein concentration by RIA and/or Lowry

Storage / Shipping

Product is shipped using dry ice (*not temperature monitored*). On arrival, store frozen below -60°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.