



PRODUCT SPECIFICATION

Thyroid Stimulating Hormone, α -Subunit (α -TSH), Human

Catalog Number: T0214 **Part Number:** 90381

Source: Human Pituitary Glands

Form

White to off-white powder

Purity Specification

≥ 90% by SDS-PAGE

Report contaminant levels if available

Characteristics Specification

This material has been bulk processed and is bottled by dry weight.

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. This product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.