



PRODUCT SPECIFICATION

Troponin C, Human Heart

Catalog Number: T5224 **Part Number:** 90338

Source: Human Heart

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.02% Sodium Azide, pH 7.4 ± 0.1

Purity Specification

≥ 95% by SDS-PAGE

Protein / Content Specification

≥ 0.5 mg/ml by the Bradford method (BioRad Protein Assay, BSA standard)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.