



Certificate of Conformance

MATERIAL

| | |
|---------------------|---|
| Description | Buffer, Conjugate Dilution |
| Catalog No. | B0221 |
| Part No. | 90092 |
| Lot No. | 2260001 |
| Storage | 2°C to 8°C. Do Not Freeze! |
| Form | Clear yellow-golden liquid; buffer containing proprietary stabilizing ingredients. |
| Intended Use | For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures. |
| Precaution | This product is not considered a biohazard. Handle using normal laboratory procedures. The buffer should not contact eyes, skin, or clothing. Wear gloves, lab coat, and eye protection when handling. Wash hands thoroughly after use. |

BATCH DATA

General Information:

This product is designed for the resuspension and dilution of colloidal gold conjugates. The buffer promotes the stabilization of colloidal gold conjugates dried on the conjugate matrix, and provides for rapid reconstitution.

Expiration Date: August 7, 2008

Manufacturing/Assay Date - 2/7/2007

On-line Certificate

Signature Not Required May 21, 2007

Laboratory Manager

PURIFIED BIOCHEMICALS AND ANTIBODIES • IMMUNOASSAY DEVELOPMENT

6838 Flanders Drive, San Diego, CA 92121-2904 USA Tel: 858-546-5800 Fax: 858-546-5812 E-Mail: contact@scrippslabs.com www.scrippslabs.com