



Certificate of Conformance

Material IgE, Purified (myeloma, kappa & lambda)
Catalog Number I0224
Scripps Part Number 90140
Lot Number 2131302
Source Human Serum (Myeloma, kappa & lambda)
Form Clear liquid - Buffer: 10mM Potassium Phosphate, 100mM Sodium Chloride, 0.1% Sodium Azide, pH 8.0 ± 0.1
Storage 2°C to 8°C short term (up to 2 weeks); below -15°C long term. Avoid freeze-thaw cycles.
Intended Use For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.
Precaution **Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.**
Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Product Data

Specifications	Results
Purity	
≥ 99% by SDS-PAGE	> 99%
Protein /Content	
≥ 1mg/ml by Lowry Protein Assay (BSA standard)	3.5 mg/ml
Activity	
Report activity (IU/ml) and method	1,120,000 IU/ml by ACS:180
Bayer ACS:180 is calibrated against WHO 2 nd IRP 75/502	

Q.C. Verification

*On-Line Certificate
Signature Not Required* 8/22/2003

Shipping Verification

Laboratory Manager

Order Processing Associate