



# Certificate of Conformance

**Material** IgE Plasma  
**Catalog Number** I0323  
**Scripps Part Number** 90141  
**Lot Number** BA247002  
**Source** Human Plasma  
**Form** Liquid, neat plasma without preservative  
**Storage** 2°C to 8°C short term (up to 2 weeks); below -15°C long term. Avoid freeze-thaw cycles.  
**Intended Use** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.  
**Precaution** **Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.**  
Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

## Product Data

Specifications	Results
	<b>Purity</b>
Neat human plasma	Neat human plasma with no preservative
	<b>Protein /Content</b>
Report IgE concentration (2.4ng/IU)	> 36 µg/ml
	<b>Activity</b>
Report activity (IU/ml)	> 15,000 IU/ml
Activity is based on WHO 2 <sup>nd</sup> IRP 75/502, where 1 IU = 2.4 ng IgE	

### Q.C. Verification

*On-Line Certificate* 12/22/2003  
*Signature Not Required*

### Shipping Verification

---

Laboratory Manager

---

Order Processing Associate