



Certificate of Conformance

Material Immunoglobulin M (IgM), purified
Catalog Number I1124
Scripps Part Number 90145
Lot Number 2171901
Source Human Plasma
Form Clear liquid - Buffer: 50mM Tris, 200mM Sodium Chloride, 0.1% Sodium Azide, pH 8.0 ± 0.1
Storage 2°C to 8°C short term (up to 1 week); below -15°C long term. Avoid freeze-thaw cycles.
Intended Use For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.
Precaution **Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.**
Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Product Data

Specifications	Results
Purity	
Report purity by immunoelectrophoresis (IEP)	Single arc by IEP using antibody to whole human serum
Protein /Content	
Report protein concentration by A_{280} ($E_{mg/ml} = 1.33$)	1.48 mg/ml

Q.C. Verification

*On-Line Certificate
Signature Not Required* 10/12/2004

Shipping Verification

Laboratory Manager

Order Processing Associate