



Certificate of Conformance

Material Immunoglobulin G (IgG), purified
Catalog Number I1424
Scripps Part Number 90146
Lot Number BA188007
Source Human Plasma
Form Clear liquid - Buffer: 20mM Sodium Phosphate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.4 ± 0.1
Storage 2°C to 8°C short term (up to 1 week); below -15°C long term. Avoid freeze-thaw cycles.
Intended Use For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.
Precaution **Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.**
Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Product Data

Specifications	Results
	Purity
Report purity by immunoelectrophoresis (IEP)	Single precipitin arc using antisera to whole human serum. No cross reaction with antisera to human IgA, IgD, IgE, or IgM
	Protein /Content
Report protein concentration by A_{280} ($E_{mg/ml} = 1.35$)	24.3 mg/ml

Q.C. Verification

On-Line Certificate 7/17/2003
Signature Not Required

Shipping Verification

Laboratory Manager

Order Processing Associate