



Certificate of Conformance

Material Immunoglobulin A (IgA), purified
Catalog Number I2224
Scripps Part Number 90147
Lot Number BF283002
Source Human Plasma
Form Clear liquid - Buffer: 20mM Tris, 100mM Sodium Chloride, 0.1% Sodium Azide, pH 8.0 ± 0.1
Storage 2°C to 8°C short term (up to 1 week); below -15°C long term. Avoid freeze-thaw cycles.
Intended Use For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.
Precaution **Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.**
Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Product Data

Specifications	Results
	Purity
Report purity by immunoelectrophoresis (IEP)	Single precipitin arc using antisera to whole human serum. No cross reaction with antisera to human IgD, IgE, IgG, or IgM
	Protein /Content
Report protein concentration by A_{280} ($E_{mg/ml} = 1.35$)	2.23 mg/ml

Q.C. Verification

*On-Line Certificate
Signature Not Required* 7/17/2003

Shipping Verification

Laboratory Manager

Order Processing Associate