



Certificate of Conformance

Material Immunoglobulin G, Fc Fragment (IgG, Fc), purified
Catalog Number I2724
Scripps Part Number 90148
Lot Number 2143801
Source Human Plasma
Form Clear liquid - Buffer: 10mM Sodium Phosphate, 250mM Sodium Chloride, pH 7.6 ± 0.1. Material has been sterile filtered and bottled in autoclaved bottles. Product does not contain any preservative; handle aseptically.
Storage Below -15°C. Freeze immediately after aliquoting.
Intended Use For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.
Precaution **Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.**
Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Product Data

Specifications	Results
Purity	
Report purity by immunoelectrophoresis (IEP)	Single arc by IEP using anti-whole human serum and anti-human IgG, Fc. No reaction to anti-IgG, Fab
Protein /Content	
Report protein concentration by A ₂₈₀ (E _{mg/ml} = 1.4)	9.8 mg/ml by A ₂₈₀

Q.C. Verification

*On-Line Certificate
Signature Not Required* 12/17/2003

Shipping Verification

Laboratory Manager

Order Processing Associate