



## Certificate of Conformance

**Material** IgE Myeloma Serum, Lambda  
**Catalog Number** I3323  
**Scripps Part Number** 90150  
**Lot Number** 2162702  
**Source** Human Serum (Myeloma)  
**Form** Liquid, neat serum  
**Storage** 2°C to 8°C short term (up to 2 weeks); below -15°C long term. Avoid freeze-thaw cycles.  
**Intended Use** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.  
**Precaution** **Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.**  
Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

### Product Data

Specifications	Results
	<b>Purity</b>
Neat human serum	Neat human serum
	<b>Protein /Content</b>
Report IgE concentration	> 40 mg/ml
	<b>Activity</b>
Report activity (IU/ml) and method	> 18,000,000 IU/ml

#### Q.C. Verification

*On-Line Certificate  
Signature Not Required* 6/14/2004

#### Shipping Verification

---

Laboratory Manager

---

Order Processing Associate