

Certificate of Conformance

MATERIAL

Description	Goat antibody to Human Chorionic Gonadotropin, α -Subunit (α -hCG), purified
Catalog No.	GC089
Part No.	90121
Lot No.	2150201
Source	Goat Serum
Immunogen	Purified Human Chorionic Gonadotropin, α -Subunit
Purified By	Affinity chromatography on immobilized hCG
Storage	2°C to 8°C. Product may be frozen below -15°C. Avoid freeze-thaw cycles.
Form	Clear liquid, Buffer: 10mM Sodium Phosphate, 500mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 \pm 0.1
Intended Use	For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.
Precaution	This product is not considered a biohazard. Handle using normal laboratory procedures.

BATCH DATA

Purity Specification:

Affinity purified, free from particulates.
Report purity by SDS-PAGE

Protein/Content Specification:

≥ 5 mg/ml by A_{280} ($E_{mg/ml} = 1.38$)

Purity Result:

> 95% by SDS-PAGE

Protein/Content Result:

7.93 mg/ml

Manufacturing/Assay Date - 2/19/2004

On-line Certificate

Signature Not Required May 31, 2007

Laboratory Manager